

# **SLIDE-IN CASSETTE FOR A CUP FOR TESTING OF DRUGS OF ABUSE**

## **1. FIELD OF THE INVENTION**

This invention relates to the art of handling, testing, and transporting fluid specimens. More particularly, it relates to a cup with a slide-in cassette to provide testing of drugs of abuse in bodily fluids, such as urine, blood, saliva, etc.

## **2. BACKGROUND OF THE INVENTION**

Fluid specimens, particularly urine, are normally collected in containers, vials or cups. When it is desired to run tests on liquid or fluid specimens contained in the cups, the lids are normally removed and specimen samples are taken out of the cups and transferred to a test apparatus. In the Instacheck® Drug Screen Drug Test, a urine sample from a cup is drawn up in a pipette and 3-4 drops (~0.2 ml) are then dispensed onto the sample well. The urine then travels up a chemical strip for 3-8 minutes. The chemical strip was pre-coated with drug conjugate on the test band. A colored anti-drug monoclonal antibody colloidal gold conjugate pad is placed at one end of the strip. In the absence of the drug in the urine, the colored antibody colloidal gold conjugate moves along the sample solution upward on the strip chromatographically by the capillary action to the immobilized drug conjugate zone on the test band region and attaches to the drug conjugate to form a visible line on the antibody complexes with the drug conjugate. Therefore, the formation of a visible precipitate in the test zone occurs when the test urine is negative for the drug. When drug is present in the urine, the drug/metabolite antigen competes with drug conjugate on the test band region for the limited antibody sites on the antibody-colloidal gold conjugate. When a sufficient concentration of drug is present, it will fill the limited antibody binding sites. This will prevent attachment of the colored antibody-colloidal gold conjugate to the drug conjugate zone on the test band region. Therefore, absence of the color band on the test region indicates a positive result.

A difficulty with the Instacheck® test is that the urine needs to be transferred from a cup onto test strips with the lid of the cup removed, thus exposing the operator and work area to possible contamination.

Additionally, the specimen sample could become contaminated as well as the worker and the surrounding equipment. Furthermore, with lid removed, spillage and loss of the unique specimens may occur. Thus, it is the object of this invention to provide a custom designed integrated system composed of a custom collection cup used as a collection and testing vessel and a custom designed slide-in test cartridge to test for drugs of abuse and other chemical and biological substance in urine and other liquid mediums in a closed, safe and secure environment.

## 2. Description of the Prior Art

U. S. Patent No. 5,119,830 to Davis describes a specimen cup having a valve to selectively operated from outside the specimen cup to introduce fluid specimen for detection of drugs of abuse by chemical strips.

U. S. Patent No. 5,916,815 to Lappe describes a specimen cup to detect drugs of abuse using intentional false positive to initially preserve anonymity.

U. S. Patent No. Des. 404,812 describes a multiple drug test card to be housed in a cup for detection of drugs of abuse. It requires sliding a card through a slotted lid and thus exposure, spillage and contamination are possible. The card is neither sealed nor contained within the device and thus can contaminate specimen. Additionally, the card draws sample from the side and required both a maximum and minimum fill requirement which makes exposure and spillage a greater problem as user tries to fill container "just right". If the minimum and maximum fill marks are not followed the test will not function. Too little urine and the test does not run, too much and the test sample is contaminated. The card must be removed at the completion of the test cycle, resulting in exposure and contamination to user and work area. If the sample is positive, the cover is removed and a closed cover is placed on bottle. Again, exposure and spillage is a problem. Lastly again, the card is inserted in the middle of a low bottle resulting in difficulty

in reading result and often requiring the user to lift the card out to view or tip bottle to view. Either way exposure and spillage is a problem.

All of the above patents had to use very complicated and/or expensive collection/reagent system. They are troublesome to get quick and easy test results. Additionally, some result in difficulty in transporting or storing the fluid specimen.

## **SUMMARY OF THE INVENTION**

Accordingly, it is an object of the present invention to provide a easy to use, inexpensive, integrated testing system comprised of a collection cup/testing vessel and a slide-in testing cassette housing the chemical/immunological test strips for the testing of drugs of abuse and other chemical and biological substances in urine and other liquid specimens/samples. The integrated system is composed of the custom test cup used to collect the sample and then the same cup is used as the testing vessel and ultimately as the storage and transport container. The test cup also can comprise a spill-prevention and over-fill prevention flap or float. This component is a movable device that is in a vertical position at the start of filling. As urine or other liquid sample is placed in the cup the "flap" will raise to a horizontal position. When raised it cuts off the available space in the cup and creates an artificially filled environment preventing additional liquid from being added to the cup. The cup is designed with a "flat" face, set back in the circular cup to move the viewing area closer to the test device while maintaining a circular type cup at the top and bottom for stability and ease of use. The "flat" viewing window also results in a ergonomically designed cup that is easier to handle when the subject is providing the urine or other sample. The inside bottom of the cup is designed with a sloped bottom (1-3 degrees) to allow for the urine sample or other liquid sample to be channeled towards the test cassette, thus allowing for testing when small volumes of specimen are given. The test cassette is uniquely designed to draw urine from the bottom, thus minimizing the amount of urine needed to perform the test. This design also eliminates the need for minimum sample volume requirements or having to tilt, turn or

invert the container to allow sample to contact the test strips. The card is hermetically sealed both around the entire perimeter as well as vertically between each test strips and horizontally below the test regions. This assures that each test strip is isolated within a unique test column and prevents any cross-contamination between the chemicals/substances contained within each test strip. The area of the card where the test regions of the test strips are viewed is covered with a clear material hermetically sealed to the face of the test card to prevent any direct contamination of the test strips from the sample or tampering with the test strips by the operator or donor. There is a sample "pooling" area at the bottom of the test cassette to allow urine or other liquid sample to migrate up to contact the test strips. This "pooling" area functions as an internal sample well. This allows the test strips to be completely enclosed in the device and eliminates any contact from the operator or donor which could cause contamination. Additionally, running horizontally above the "pooling" area is a "dam" designed to restrict the vertical flow of sample up the test strips and contain the sample in the "pooling" area.

#### **BRIEF DESCRIPTION OF THE DRAWINGS**

The novel features which are believed to be characteristic of the invention are set forth with particularity in the appended claims. The invention itself, however, both as to its organization and method of operation, may best be understood by reference to the following description, when taken in connection with the accompanying drawing in which:

FIG.1 is a prospective view of the cup for testing drugs of abuse and other chemical and biological substance of the present invention.

FIG.2 is a top prospective view of the specimen cup (with the top cover removed) of FIG.1; and

FIG.3 is a (bottom) prospective view of the slide-in cassette.

## **DETAILED DESCRIPTION OF THE PREFERRED EMBODIMENT**

5 A specimen cup (100) of the present invention includes a base container (112) and a lid (114). The specimen cup (100) is for collecting, testing, storing and transporting a urine specimen and other liquid within a container thereof. The base container (112) can optionally has an expanded sample collection portion to allow more urine to be collected. The lid (114) has threads (116) which mesh with threads (118) of the base container (112) to sealingly hold the lid (114) on the base container 112. In region A behind the chamber (104) to which the cassette (102) is hermetically sealed, there is a urine spill prevention flap or float (108) (see FIG. 2) to which the urine once entered into the sample collection portion will be prevented from  
10 splashing during transport or storage. The flap or float is free to travel vertically in region A under the pressure from the fluid specimen, such as urine.

15 The base container (112) and its lid (114) are constructed of a material which is transparent, and impervious to fluid specimens contained therein. The materials include but not limited to thermoplastics, specialty plastics, thermosets, engineering plastics.

20 Thermoplastics include but not limited to: polyamideimide (PAI), polyethersulfone (PES), polyarylsulfone (PAS), polyetherimide (PEI), polyarylate (PAR), polysulfone (PSO), polyamide (PA), polycarbonate (PC), styrene-maleic anhydride (SMA), chlorinated PVC (CPVC), poly(methylmethacrylate) (PMMA), styrene-acrylonitrile (SAN), polystyrene (PS), acrylonitrile-butadiene-styrene (PS), acrylonitrile-butadiene-styrene (ABS), poly(ethyleneterephthalate) (PET),  
25 poly(vinylchloride) (PVC), polyetherketone (PEK), polyetheretherketone (PEEK), polytetrafluoroethylene (PTFE), poly(phenylene sulfide) (PPS), liquid crystal polymer (CCP), nylon-6,6, nylon-6, nylon-6,12, nylon-11, nylon 12, acetal resin, low and high density polypropylene (PP), high  
30 density polyethylene (HDPE), low density polyethylene (LDPE), polystyrene, ethylene-vinyl acetate, poly-vinyl-acetate, polyacrylic, etc., or a copolymer or a combination thereof.

Specialty plastics include but not limited to fluorocarbon polymers  
 5 and infusible film products such as Kapton, Upilex polyimide film etc., a  
 copolymer or a combination thereof. Thermosets include but not limited to  
 phenolics, epoxies, urea-formaldehyde, silicones, etc., a copolymer or a  
 combination thereof. Engineering plastics include but not limited to acetyl  
 10 resins, polyamide, polyetherimides, polyesters, liquid crystal polymers,  
 polycarbonate resins, poly(phenylene ether) alloys, polysulfone resins,  
 polyamideimide resins, etc., a copolymer or a combination thereof .

The bottom floor (120) of the cup can be optionally sloped from the  
 backside (122) downwardly at 1-3° towards the front side (124). This  
 15 forces the fluid (by gravity) to moves forward, hence reduces the fluid  
 specimen needed for the testing for drugs of abuse by the cassette. The  
 front of the cup has a retracted flat face (200) designed to move the viewing  
 area closer to the test cassette. The base and top of the cup remain circular  
 to allow for use of standard covers and provide a stable base. Inside the cup  
 20 are custom channels (156) used to guide and oriented the cassette in the  
 device. The cassette is inserted into the cup (100) with its outside edges  
 (150) anchored between the bars (158). The slot on the left side of the  
 cassette will only align with the triple channel on the left side of the cup.  
 The bars (158) ensure the cassette is inserted facing the correct way for  
 25 viewing and ensure proper placement within the container. Because one of  
 the fluids that may be tested is urine, as the urine cools in a closed  
 environment condensation may occur. The tracks are designed to orient the  
 cassette for viewing while allowing movement of air between the cassette  
 and face of cup to prevent condensation forming on inside of cup. The  
 30 chemical test strips (106) of various, flexible configurations such as 11-nor-  
 $\Delta$ -9-tetra hydrocannabinol-9 carboxylic acid (THC), Cocaine (COC),  
 Methamphetamine/amphetamine (MAP), 1-(1'- phenylcyclohexyl)  
 piperidine (PCP), Morphine (MOR) etc. are housed in a custom cassette  
 (126). The cassette has four distinct, isolated test channels (132, 134, 136  
 35 and 138) which house the test strips. Each test channel has a clear, sealed  
 window for viewing the results. Each channel is hermetically sealed both  
 vertically and horizontally

to ensure four unique test areas and prevent any direct or cross  
 5 contamination. As seen in FIG. 3 the cassette is formed by an upper (128)  
 and lower (130) member. Near the bottom of the cassette is a horizontally  
 running "dam" (260) that when the upper and lower members are  
 hermetically sealed together creates a sample "pooling" area (210). This  
 "pooling" area (210) allows sample to contact the test strips while  
 10 eliminating the need for the test strips (106) to be exposed. Thus the entire  
 test strip is contained within the cassette eliminating potential  
 contamination, adulteration or tampering. When the test card is inserted  
 into the test container, the sample "pools" around the base of the test strips  
 and wicks vertically up the strips. As the sample moves up the strip, the  
 15 result is observed through the clear viewing windows. The clear viewing  
 windows prevent direct contact with the test regions of the test strips either  
 by the operator, donor or specimen.

During operation, a specimen, such as urine, is provided in the  
 20 custom collection/test cup (100). The test cassette (126) is inserted into the  
 test cup through custom bars (156) and the lid of the cup (114) is put in  
 place. The urine specimen then enters the "pooling" areas (210) at the base  
 of the test cassette and begins to wick up the test strips. When the urine  
 contacts the test strips, the characteristics thereof, in conjunction with  
 25 chemicals in the test strips causes the test strips to change color, thereby  
 providing a visual indication to an operator in accordance with the  
 precalibrated indicator marking beside the respective test strips  
 corresponding to such characteristics. The changes in color are then easily  
 observed and read by the operator through the transparent window on the  
 30 test card and the face of the collection/test cup.

After testing is completed, the specimen can be stored, transported or  
 disposed of in the collection/test cup used for this testing process. This  
 eliminates having to remove the test device, change lids, transfer specimen  
 35 or otherwise handle the urine sample in any way that could result in  
 exposure or contamination to the operator, donor or surrounding  
 environment.